

Ultrasensitive, Multiplexed Detection of Cancer Biomarkers Directly in Serum by Using a Quantum Dot-Based Microfluidic Protein Chip

Mei Hu, Juan Yan, Yao He, Haoting Lu, Lixing Weng, Shiping Song, Chunhai Fan Lianhui Wang

Laboratory of Advanced Materials, Fudan University Laboratory of Physical Biology, Shanghai Institute of Applied Physics, Functional Nano and Soft Materials Laboratory, Suzhou University,



ACS NANO ASAP



Introduction


1. Semiconductor quantum dots (QDs) and its fluorescent tuning have attracted intensive attention during the last two decades.
2. The QDs, which are synthesized in organic phase, often possess a high photoluminescent quantum yield (PLQY 60-85%) and narrow size distribution.
3. QDs that are directly prepared in the aqueous phase (aqQDs) possess excellent aqueous dispersibility because their stabilizers (*e.g.*, thiols) are naturally water-dispersible and can be used for biological studies directly.
4. The QD-based bioconjugates were already applied into ELISA, Western blotting, and microarray for protein assays.
5. However, all of these techniques are based on a static solid/liquid interface reaction. In such systems, QDs are less active than a small molecular dye because of steric hindrance and are prone to deposit on surfaces.
6. To address this challenge, microfluidic chips based on the manipulation of a continuous liquid flow through microfabricated channels are introduced in this work.
7. Herein this group report the synthesis of antibody-conjugated aqQDs with strong photoluminescence and robust stability and an aqQD-based microfluidic protein chip for ultrahigh sensitive, selective, and multiplex detection of cancer biomarkers.


Preparation of QD-IgG Conjugates (anti-mouse igG)

First, waterdispersed QDs (in 150 mM PBS, pH 7.3) were activated with EDC/NHS. To 60 L of QDs (12 M), add 17.25 L of EDC (33.4 mM in H₂O) and 16.25 L of NHS (70.9 mM in H₂O). The reaction solution was incubated at 25 °C for 15 min with gentle shaking. EDC is easily hydrolytic and should be dissolved in water immediately before use. Then 27 L of goat anti-mouse IgG (13.3 M) in PBS buffer was added to the activated QDs. The molar ratio of QD/EDC/NHS/IgG was calculated as 1:800:1600:0.5.



a**b**





Summary

- 1. The present study has successfully conjugated QDs with a secondary antibody.**
- 2. This bionanohybrid exhibits excellent performance as a fluorescent probe in protein assays.**
- 3. The QD-IgG conjugate, combined with a microfluidic protein chip, improved the detection limits of cancer biomarkers up to 250 fM, which represents a sensitivity improvement up to 4 orders of magnitude as compared to organic dyes (*e.g.*, *commercially* available Cy3-Abs conjugates).**
- 4. Besides, multicolored imaging was realized and potentially useful for multiplexed detection in proteomic studies. This study thus opens a new path for practical clinical and proteomic applications of protein chips.**

Thanks