

# Immunohistochemistry for the detection of neural and inflammatory cells in equine brain tissue

Gretchen H. Delcambre<sup>1,\*</sup>, Junjie Liu<sup>2</sup>, Jenna M. Herrington<sup>2</sup>, Kelsey Vallario<sup>2</sup> and Maureen T. Long<sup>2,\*</sup>

<sup>1</sup> Department of Biomedical Sciences/College of Veterinary Medicine and Biomedical Sciences, Colorado State University, Fort Collins, Colorado, USA

<sup>2</sup> Department of Infectious Diseases and Pathology/College of Veterinary Medicine, University of Florida, Gainesville, Florida, USA

\* These authors contributed equally to this work.

## ABSTRACT

Phenotypic characterization of cellular responses in equine infectious encephalitis has had limited description of both peripheral and resident cell populations in central nervous system (CNS) tissues due to limited species-specific reagents that react with formalin-fixed, paraffin embedded tissue (FFPE). This study identified a set of antibodies for investigating the immunopathology of infectious CNS diseases in horses. Multiple commercially available staining reagents and antibodies derived from antigens of various species for manual immunohistochemistry (IHC) were screened. Several techniques and reagents for heat-induced antigen retrieval, non-specific protein blocking, endogenous peroxidase blocking, and visualization-detection systems were tested during IHC protocol development. Boiling of slides in a low pH, citrate-based buffer solution in a double-boiler system was most consistent for epitope retrieval. Pressure-cooking, microwaving, high pH buffers, and proteinase K solutions often resulted in tissue disruption or no reactivity. Optimal blocking reagents and concentrations of each working antibody were determined. Ultimately, a set of monoclonal (mAb) and polyclonal antibodies (pAb) were identified for CD3<sup>+</sup> (pAb A0452, Dako) T-lymphocytes, CD79α<sup>+</sup> B-lymphocytes (mAb HM57, Dako), macrophages (mAb MAC387, Leica), NF-H<sup>+</sup> neurons (mAb NAP4, EnCor Biotechnology), microglia/macrophage (pAb Iba-1, Wako), and GFAP<sup>+</sup> astrocytes (mAb 5C10, EnCor Biotechnology). In paraffin embedded tissues, mAbs and pAbs derived from human and swine antigens were very successful at binding equine tissue targets. Individual, optimized protocols are provided for each positively reactive antibody for analyzing equine neuroinflammatory disease histopathology.

Submitted 2 October 2015  
Accepted 23 December 2015  
Published 25 January 2016

Corresponding author  
Gretchen H. Delcambre,  
gh.delcambre@colostate.edu

Academic editor  
María Ángeles Esteban

Additional Information and  
Declarations can be found on  
page 12

DOI 10.7717/peerj.1601

© Copyright  
2016 Delcambre et al.

Distributed under  
Creative Commons CC-BY 4.0

**OPEN ACCESS**

**Subjects** Veterinary medicine, Infectious diseases, Pathology, Histology

**Keywords** Immunohistochemistry, Equine, Neuropathology, Leukocytes, Neuroglia

## INTRODUCTION

Comparative pathology of many diseases that affect multiple species is hindered by the lack of species-specific reagents. Although many antibodies that react with cellular antigens and inflammatory molecules are available, there still are gaps in how these work in formalin-fixed paraffin embedded (FFPE) tissues. This limits our ability to generate

new knowledge about disease processes that can benefit human and animal health as a whole. Horses develop many of the same or similar central nervous system (CNS) infections as humans (*Ritchey et al., 2006; Bourgeois et al., 2011; Rushton et al., 2013; Yu et al., 2015*). Cell marker panels are often composed of both equine and non-equine specific antibodies, of which most are used in flow cytometry. Since manufacturers may not have supporting technical documentation on whether their products will cross-react with equine antigens (*Beckstead, 1994; Ramos-Vara, 2005*), development of antibody panels to accomplish in situ disease characterization in formalin-fixed tissue is a formidable task because cross-linking of antigens often renders epitopes non-reactive. Additional investigative steps are required to determine if antigens can be retrieved and what retrieval method works best. Because of this and other hurdles, there is limited information on what antibodies react with FFPE tissue.

The goal of this study was to identify a panel of cell markers for studying cellular pathogenesis in equine infectious brain diseases. In particular, protocols for phenotyping reactive glia, neurons, and infiltrating peripheral blood mononuclear and polymorphonuclear cells in the equine brain were developed. Commercially available antibodies for CD3<sup>+</sup> T-lymphocytes, CD8<sup>+</sup> and CD4<sup>+</sup> T cell subpopulations, B lymphocytes, macrophages, microglia, astrocytes, and neurons were investigated for reactivity in normal and diseased FFPE tissue. Optimization of manual immunohistochemistry (IHC) protocols was performed using various IHC staining reagents and methods.

## MATERIALS AND METHODS

### Tissue samples

Immunohistochemistry protocols were developed on diseased and normal horse tissues. Diseased horses consisted of clinically affected West Nile virus (WNV), both naturally and experimentally infected, and *Sarcocystis neurona* infections (*Beckstead, 1994; Gutierrez et al., 1999; Porter et al., 2003; Seino et al., 2007*). Neurologically, normal horses were obtained by owner surrender for humane euthanasia due to loss of use. Brain, spinal cord, lymph node, spleen, thymus, and liver were collected and archived from these animals under University of Florida (Gainesville, FL) Institutional Animal Care and Use Committee protocols #F077, #F093, #D163, and #4109. Tissues were fixed in 10% neutral buffered formalin and processed into paraffin-embedded blocks approximately one-week after fixation. Initial evaluation of antibody binding was tested on non-infected equine lymph node and spleen for lymphocytic targets, liver and thymus for tissue macrophage targets, and brain for microglia, astrocytes, and neurons.

### Tissue processing

The invariable IHC procedures for all protocols included sectioning FFPE tissues at 5  $\mu\text{m}$  and placing them on positively charged glass slides. The slides were soaked in xylene (Fischer Scientific, Pittsburg, PA, USA) three times for 5 min to remove paraffin. These sections were then rehydrated through a gradient of ethanol (Fischer Scientific) for 5 min in each concentration, 100%, 100%, 95%, and 70% ethanol, followed by de-ionized

water. In order to reduce the volume of the reagents tested and liquid loss, tissues were encircled with a hydrophobic barrier pen (ImmEdge™ Pen, Ted Pella Inc., Redding, CA, USA).

### Antigen unmasking

Three methods of heating slides for investigating heat induced epitope retrieval (HIER) effectiveness included using a pressure cooker, microwave, and a double boiler. Pressure cooking was performed at 125 °C for 30 s followed by 90 °C for 10 s (*Matyjaszek et al., 2009*; *Grosche et al., 2012*), or microwaving was performed for 10 min (*Kumar & Rudbeck, 2009*) using a commercial countertop GE 1000W oven. For double boiling, two tissue-slides were floated back-to-back in 25 ml of retrieval solution in a 50 ml plastic conical tube. Conical tubes were then placed in pre-warmed water of a 250 ml glass beaker on a hotplate. Water temperature was maintained at 90 °C. After approximately 5 min of warming the retrieval solution, HIER was timed for 10 min. Conical tubes were then removed from the double boiler and allowed to cool for 15 min at 27 °C. Tissues were rinsed in deionized water three times for 2 min.

Heat induced epitope retrieval buffers were tested primarily using the double boiler system. Regents included two commercial citrate buffers, Epitope Retrieval Solution pH 6 (Novocastra, Leica, Newcastle Upon Tyne, UK) and Target Retrieval Solution pH 6 (Dako, Glostrup, Denmark), and an ethylenediaminetetraacetic acid (EDTA) solution buffered at pH 9 (10 mM Tris Base, 1 mM EDTA solution, 0.05% Tween 20, and NaOH to titrate to pH 9). A 1 × concentration of each solution was freshly made by diluting stock solutions with deionized water. For proteolytic epitope retrieval, tissues were treated with 200 ug/ml proteinase K solution (Tris HCL 100 mM pH 8.2, Tween 20, and Proteinase K (Ambion, Foster City, CA, USA)) for 10 min at 37 °C.

### Endogenous peroxidase blocking

Peroxidase neutralizing solutions that were tested included two prepared solutions of hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>), 3% and 0.3% H<sub>2</sub>O<sub>2</sub>, and a ready-to-use commercial reagent, Peroxidase Block (Novolink™ Polymer Detection System; Leica, Wetzlar, Germany). Solutions containing 3% and 0.3% H<sub>2</sub>O<sub>2</sub> were made fresh for each staining attempt by diluting 30% H<sub>2</sub>O<sub>2</sub> (Fischer Scientific) in 1 × phosphate buffered saline (PBS) (10 × PBS, Fischer Scientific). Tissues were immersed in peroxidase blocking solution for 5 min followed by two, 5 min rinses in PBS.

### Non-specific protein blocking

Non-specific blocking techniques included four commercial reagents and one lab prepared solution. Commercial reagents included 10% Normal Goat Serum (Invitrogen, Frederick, MD, USA), Protein Block (Novolink™ Polymer Detection Systems; Leica, Wetzlar, Germany), Novocastra™ IHC/ISH Super Blocking Solution (Leica), and Novocastra™ Liquid Serum, Normal Goat Serum Blocking Reagent (Leica, Wetzlar, Germany). Additionally, a 5% goat serum solution was prepared by diluting Immunopure® Goat Serum (ThermoFischer Scientific, Waltham, MA, USA) in 1 × PBS.

**Table 1** Antibodies tested for immunohistochemical reactivity in equine tissues.

Antibody	Antigen species	Host species	Isotype	Name or clone	Protein target	Source
CD3	Human	Mouse	IgG1	F7.2.38	Thymocytes, T lymphocytes, natural killer cells	Dako
CD3	Human	Rabbit	–	A0452	T lymphocytes	Dako
CD4	Equine	Mouse	IgG1	CVS4	T helpers	AbD
CD4	Human	Mouse	IgG1	4B12	Thymocytes, T helpers, mature peripheral T cells, TReg or cytotoxic T cell subsets	Dako
CD4	Human	Mouse	IgG1	1F6	T helpers	Leica
CD4	Equine	Mouse	IgG1	HB61A	Thymocytes, T helpers	VMRD
CD8	Equine	Mouse	IgG1	CVS8	Cytotoxic T lymphocytes	AbD
CD8	Equine	Mouse	IgG1	HT14A	Thymocytes, Cytotoxic T lymphocytes	VMRD
CD5	Equine/Bovine	Mouse	IgG2A	B29A	B lymphocytes only, in horses	VMRD
CD20	Human	Mouse	IgG1	MJ1	B lymphocytes	Leica
CD21	Human	Mouse	IgG1	Bu33	B lymphocytes	AbD
CD21	Human	Mouse	IgG1	1F8	B lymphocytes	Dako
CD79 $\alpha$ cy	Human	Mouse	IgG1	HM57	B lymphocytes	Dako
IgG (H&L)	Equine	Goat	IgG	A70–106	B lymphocytes	Bethyl
IgG (H&L)	Equine	Rabbit	IgG	A70–118	B lymphocytes	Bethyl
IgG (H&L)	Goat	Rabbit	IgG	A50–100	B lymphocytes	Zymed
Macrophage	Rabbit	Mouse	IgG1	RAM11	Macrophages	Dako
Macrophage	Human	Mouse	IgG1	MAC387	Macrophages	Leica
Macrophage	Human	Mouse	IgG1	AM-3K	Macrophages	TransGenic
CD68	Human	Mouse	IgG1	CD68-KP1	Macrophages	Leica
NF-H	Bovine	Mouse	IgG2b	9B12	Neurofilament-Heavy chain	EnCor-Biotechnology
NF-H	Bovine	Mouse	IgG1	AH1	Neurofilament-Heavy chain	EnCor-Biotechnology
NF-H	Swine	Mouse	IgG1	NAP4	Neurofilament-Heavy chain	EnCor-Biotechnology
GFAP	Swine	Mouse	IgG1	5C10	Astrocytes, Satellite cells, Schwann cells	EnCor-Biotechnology
Microglia	Human	Rabbit	–	Iba-1	Microglia, macrophages	Wako
Neu-N	Mouse	Mouse	IgG1	A60	Neurons	Millipore

Protein blocking reagents were applied for 20 min with the exception of Novolink™ Protein block for only 5 min. Blocking reagents were removed without rinsing before adding primary antibody.

### Primary antibodies

Primary antibodies tested targeted cell populations of astrocytes, microglia, neurons, T lymphocytes, B lymphocytes, and macrophages (Table 1). All primary antibodies were diluted in one of two commercial diluents (IHC Diluent (Novacastra, Leica); Dako Antibody Diluent (Dako, Glostrup, Denmark)). Two-fold serial dilutions of each antibody were tested to determine an optimal staining range. If the signal was weak to absent or background staining was present, additional dilution tests were performed until optimal staining was achieved. Antibodies were applied in a 37 °C, humidified incubator for

60–120 min or overnight (~16 hrs) at 4 °C in a humidified, covered dish. After primary antibody incubation, slides were washed three times for 5 min each in 1 × PBS before secondary antibody application.

Negative controls for each primary antibody consisted of either an isotype-matched negative primary control (MCA928, AbD Serotec, Kidlington, UK) for monoclonal antibodies or rabbit serum for polyclonal antibodies.

### Detection systems

Three commercial, horseradish peroxidase (HRP)-linked conjugate detection kits were utilized by following manufacturer's instructions. These kits included the Vectastain<sup>®</sup> ABC Kit–Mouse IgG (Vector Laboratories, Burlingame, CA, USA), the Vectastain<sup>®</sup> ABC Kit–Rabbit IgG, and the Novolink<sup>™</sup> Polymer Detection System (Leica, Wetzlar, Germany).

Lastly, a substrate-chromogen (Vector NovaRED Peroxidase Substrate, Vector Laboratories) was applied for 10 min, rinsed with de-ionized water for 5 min, counterstained with Lab Vision<sup>™</sup> Mayer's Hematoxylin (ThermoFisher Scientific) for 1 min, and rinsed in running tap water for 2 min. All sections were then dehydrated through an increasing gradient of ethanol, 50%, 70%, 95%, and 100%, for 2 min each. Slides were cleared in xylene three times for 3 min each before coverslipping with Permount<sup>™</sup> mounting medium (Fischer Scientific, Hampton, NH, USA).

## RESULTS

### Antigen unmasking

Of the reagents tested, low pH citrate based solutions resulted in superior staining (Table 2). Target Retrieval Solution pH 6 by Dako was sufficient for most antibodies; however, use of Epitope Retrieval Solution pH 6 by Leica resulted in positive staining when the Dako product failed to produce staining with the CD79αγ<sup>+</sup> antibody. Proteinase K solution was tested with macrophage antibodies and resulted in increased staining intensity of MAC387<sup>+</sup> antibody.

Various issues were encountered with the method of applying heat. Utilizing, one antibody, CD79αγ<sup>+</sup>, only the double boiling method resulted in consistent positive staining of lymph nodes. Pressure cooking and microwaving resulted in limited to no staining or uneven staining, respectively. In addition, tissue disruption was minimal with the double boiling method. One antibody, Iba-1<sup>+</sup>, did not require heat retrieval and any HIER resulted in non-specific background staining.

### Endogenous peroxidase blocking

In FFPE tissues, 3% H<sub>2</sub>O<sub>2</sub> solution was reliable for reducing background staining without disruption of the tissue. There was no difference between commercial and lab diluted reagents.

### Non-specific protein blocking

Multiple combinations of IHC antibody diluents and non-specific protein blocking reagents were tested. No difference in staining intensity was noted between of either

**Table 2** Selected protocols for immunohistochemical visualization of resident and infiltrative cells in normal and WNV-diseased equine CNS tissues.

Antibody	ID/clone	Tissue format	Peroxidase blocker	Antigen retrieval	Blocker	Dilution	Incubation	Commercial secondary kit
CD3	A0452	FFPE	3% H <sub>2</sub> O <sub>2</sub>	Target retrieval solution pH6 <sup>5</sup>	10% normal goat serum <sup>1</sup>	1:100	60 min at 37 °C	Vectastain <sup>®</sup> ABC Kit–Rabbit <sup>3</sup>
CD4	HB61A	FFT	0.3% H <sub>2</sub> O <sub>2</sub>	None	Protein block <sup>2</sup>	1:800	120 min at 23 °C	Post primary block & Novolink polymer <sup>2</sup>
CD8	HT14A	FFT	0.3% H <sub>2</sub> O <sub>2</sub>	None	Protein block <sup>2</sup>	1:800	120 min at 23 °C	Post primary block & Novolink polymer <sup>2</sup>
CD79αCY	HM57	FFPE	3% H <sub>2</sub> O <sub>2</sub>	Epitope retrieval solution pH6 <sup>2</sup>	5% immune pure goat serum <sup>4</sup>	1:100	90 min 37 °C	Post primary block & Novolink polymer <sup>2</sup>
Macrophage	MAC387	FFPE	3% H <sub>2</sub> O <sub>2</sub>	Target retrieval solution pH6 <sup>5</sup> or proteinase K	10% normal goat serum <sup>1</sup>	1:100	60 min at 37 °C	Vectastain <sup>®</sup> ABC Kit–Mouse <sup>3</sup>
NF-H	NAP4, AH1, 9B12	FFPE	3% H <sub>2</sub> O <sub>2</sub>	Epitope retrieval solution pH6 <sup>2</sup>	Protein block <sup>2</sup>	1:1,000	60 min at 37 °C	Post primary block & Novolink polymer <sup>2</sup>
GFAP	5C10	FFPE	3% H <sub>2</sub> O <sub>2</sub>	Target retrieval solution pH6 <sup>5</sup>	5% immune pure goat serum <sup>4</sup>	1:8,000	60 min at 37 °C	Vectastain <sup>®</sup> ABC Kit–Mouse <sup>3</sup>
Iba-1	–	FFPE	3% H <sub>2</sub> O <sub>2</sub>	None	Protein block <sup>2</sup>	1:500	60 min at 37 °C	Post primary block & Novolink polymer <sup>2</sup>

**Notes:**

FFPE, Formalin-fixed, paraffin embedded; FFT, Fresh, frozen tissue.

<sup>1</sup>Invitrogen<sup>™</sup>.

<sup>2</sup>Leica<sup>®</sup>–Novocastra<sup>™</sup> Product Line.

<sup>3</sup>Vector Labs<sup>®</sup>.

<sup>4</sup>FischerScientific<sup>®</sup>.


<sup>5</sup>Dako<sup>®</sup>.

diluent used. Among protein blocking solutions, 10% Normal Goat Serum, Protein Block from the Novolink<sup>™</sup> Polymer Detection System, and 5% Immunopure<sup>®</sup> Goat Serum in PBS solution were all equal in blocking non-specific staining.

### Specificity and sensitivity of monoclonal and polyclonal antibodies

Several antibodies were tested for lymphocyte specific CD antigens including CD3<sup>+</sup>, CD4<sup>+</sup>, CD5<sup>+</sup>, CD8<sup>+</sup>, CD20<sup>+</sup>, CD21<sup>+</sup>, and CD79αy<sup>+</sup> (Table 1). Of the two CD3<sup>+</sup>, pan-T cell markers examined, rabbit anti-human polyclonal antibody (A0452, Dako) appropriately stained lymph node cortex (Fig. 1A), periarterial lymphatic sheaths of spleen, and perivascular cuffs in WNV infected brain (Fig. 1B). Four CD4<sup>+</sup> T-helper cell antibodies were tested. Only mouse, monoclonal anti-equine antibody (HB61A, VMRD, Pullman, WA, USA) at 1:25 dilution positively stained lymph node cortex; however, background staining was high when applied to brain tissues. Two CD8<sup>+</sup> cytotoxic T cell markers were investigated, but neither marker had reactivity in FFPE tissue.

Antibodies against CD5<sup>+</sup>, CD20<sup>+</sup>, CD21<sup>+</sup>, CD79αy<sup>+</sup>, and IgG (H+L) (Table 1) were tested for identification of B cell populations. A putative lymphocyte marker, CD5<sup>+</sup> (B29A, VMRD), reportedly selects for B cells in equine tissues. Although this antibody intensely stained the germinal centers of FFPE lymph nodes, cortical staining was also noted. Because of this, CD5<sup>+</sup> (B29A) was unreliable for distinction between B cell and




**Figure 1** IHC of CD3<sup>+</sup> T lymphocytes and CD79<sup>+</sup> B lymphocytes in equine tissues. For the detection of lymphocytes, (A, C) equine lymph node cortex and (B, D) WNV infected equine brain incubated with (A, B) CD3<sup>+</sup> T lymphocyte primary antibody (pAb A0452; Dako, Glostrup, Denmark) for 60 minutes at 37 °C and detected by Vectastain<sup>®</sup> ABC Kit, or incubated with (C, D) CD79 $\alpha$ cy<sup>+</sup> B lymphocyte primary antibody (mAb HM57, Dako) for 90 minutes at 37 °C and detected by Novolink<sup>™</sup> Polymer Detection System. Vector NovaRED Peroxidase Substrate chromogen and hematoxylin counterstain. Bar, 50  $\mu$ m.

T cell populations. Additionally in brain tissue, this antibody resulted in non-specific background staining, which could not be resolved. No staining was achieved with CD20<sup>+</sup>, CD21<sup>+</sup>, or any IgG (H+L) antibodies. Anti-human CD79 $\alpha$ cy<sup>+</sup> (HM57; Dako, Glostrup, Denmark) monoclonal antibody at 1:100 successfully stained lymph node germinal centers with no background staining in equine brain (Figs. 1C and 1D).

Multiple macrophage-targeting antibodies were investigated. RAM 11 (Dako, Glostrup, Denmark), AM-3K (TransGenic Inc., Kobe, Japan), and CD68<sup>+</sup> (KP1; Leica, Wetzlar, Germany) antibodies had no reactivity with control tissues; however, MAC387<sup>+</sup> (Leica, Wetzlar, Germany) was reactive. This marker positively stained control hepatic and thymic macrophages (Fig. 2A) with limited staining noted in lymph node sections. Based on cell morphology, polymorphonuclear and mononuclear cells also stained positively due to lack macrophage-specificity of the antibody. No reactivity was noted in normal equine brain. The MAC387<sup>+</sup> cell population was distinct from the distribution of CD3<sup>+</sup> lymphocytes in brain sections with inflammation due to WNV infection (Fig. 2B). Additionally, this macrophage antibody had little to no cross reactivity with brain microglia since few MAC387<sup>+</sup> cells were visualized within glial nodules of WNV<sup>+</sup> brain.

To characterize reactive gliosis, antibodies against both microglial and astrocytic markers were tested. In non-infected brain tissue, Iba-1<sup>+</sup> antibody stained scattered resident microglia with low intensity, yet microglial processes were visualized (Fig. 2C). Staining intensity was higher in WNV-infected brains (Fig. 2D). Astrocyte populations were identified with antibody against an intermediate filament, glial fibrillary acid protein



**Figure 2** IHC of MAC387<sup>+</sup> and Iba-1<sup>+</sup> macrophages and microglia in equine tissues. For the detection of macrophage/microglial lineage cells, (A) equine thymus and (B) WNV infected equine brain were incubated with MAC387<sup>+</sup> primary antibody (mAb MAC387; Leica, Wetzlar, Germany) for 60 minutes at 37 °C and detected by Vectastain<sup>®</sup> ABC Kit. Additionally, (C) normal (arrows) and (D) WNV infected horse brain were incubated with microglia/macrophage primary antibody (pAb Iba-1; Wako, Neuss, Germany) for 60 minutes at 37 °C and detected by Novolink<sup>™</sup> Polymer Detection System. Vector NovaRED Peroxidase Substrate chromogen and hematoxylin counterstain. Bar, 50 μm.


(GFAP), specific to astrocytes in CNS tissues. This antibody reacted strongly with astrocytes in both normal and infected brains with distinct astrocytic processes notable. GFAP<sup>+</sup> astrocytes were distributed within the brain parenchyma, at the glia limitans, along brain vasculature, and in areas of gliosis of WNV<sup>+</sup> brains (Fig. 3).

Three antibodies that target neurofilament heavy-chain proteins (NF-H) of neuronal axons successfully stained equine brain tissue. Swine derived, NAP4 (EnCor Biotechnology Inc., Gainesville, FL, USA) antibody, was superior and stained neurons in both normal (Fig. 4) and infected brains. Mouse monoclonal Neu-N antibody (A60, Millipore, Billerica, MA, USA), which targets neuronal perikaryon, was also tested but without successful reactivity in equine tissue.


Successful manual IHC protocols for these antibodies in equine neural tissue were identified (Table 2). It should be noted that overnight, refrigeration with CD3<sup>+</sup>, CD79<sup>+</sup>, MAC387<sup>+</sup>, GFAP<sup>+</sup>, and Iba-1<sup>+</sup> primary antibodies will successfully stain tissues and may provide flexibility in laboratory scheduling if single-day protocol completion is not possible.

## DISCUSSION

Formalin-fixed, paraffin embedded tissues are commonly archived for histological examination. Formalin fixation and paraffin embedding samples retains tissue



**Figure 3** IHC of GFAP<sup>+</sup> astrocytes in WNV infected equine brain. GFAP<sup>+</sup> astrocytes in WNV infected equine brain (A) near blood vessels and (B) at the glial limitans. IHC. Bar, 50  $\mu$ m.



**Figure 4** IHC of NF-H<sup>+</sup> neurons in equine brain. Normal equine brain incubated with NF-H<sup>+</sup> neuron primary antibody (mAb NAP4; EnCor Biotechnology, Gainesville, FL, USA), for 60 minutes at 37 °C and detected by Novolink<sup>TM</sup> Polymer Detection System. Vector NovaRED Peroxidase Substrate chromogen and hematoxylin counterstain. Bar, 50  $\mu$ m.

architecture and deactivates any potential infectious agents (*Cantile et al., 2001; Ritchey et al., 2006; van Marle et al., 2007*); however, this process can block or destroy sensitive epitopes (*Beckstead, 1994; Mérant et al., 2003; Terio et al., 2003; Ibrahim et al., 2007*). Despite reported reversal of these treatments, the tested equine antigen derived CD4<sup>+</sup> and CD8<sup>+</sup> T cell were not detectable in FFPE. This result has been described in other species (*Beckstead, 1994; Bilzer et al., 1995; Zeng et al., 1996; Gutierrez et al., 1999; Härtig et al., 2009*). Anti-CD4<sup>+</sup> (mAb HB61A, VMRD) and CD8<sup>+</sup> (mAb HT14A, VMRD) monoclonal antibodies, however, were reactive to fresh, frozen equine lymph node and infectious CNS tissue. Archiving tissues in both FFPE and fresh, frozen tissue formats is recommended when investigating epitopes that may be affected by fixation.

After successful staining of control tissues, normal and pathologic brain tissue samples were then tested alongside positive and negative tissue controls. Pathologic tissue samples included WNV-infected brain, which contains reactive gliosis and perivascular cuffing of inflammatory cells (*Cantile et al., 2001; Ibrahim et al., 2007; van Marle et al., 2007; Schnabel et al., 2013*). *S. neurona* infected tissues were also used, which contained focal nonsuppurative inflammation, mononuclear perivascular cuffing, and the occurrence of giant cells and eosinophils (*Boy, Galligan & Divers, 1990; Dubey et al., 2001*). Following this workflow of tissue testing aided in the identification of successful antibody reactivity.

Several hurdles must be overcome to optimize the interaction of antibodies with their intended targets. Aldehyde cross-links that bind tissue proteins during formalin fixation process must be removed by HIER or proteolytic epitope retrieval to allow epitopes to resume a more natural conformation and increase antibody-binding capacity (*Shi, Key & Kalra, 1991; Ferrier et al., 1998; Krenacs et al., 2010*). Peroxidases that naturally occur within tissues must be neutralized with  $H_2O_2$  in order to prevent non-specific staining during the application of peroxidase-based substrate kits (*Wendelboe & Bisgaard, 2013*). Endogenous proteins may non-specifically interact with antibodies and cause background staining that masks target antigen signal (*Daneshtalab, Doré & Smeda, 2010; Buchwalow et al., 2011*). These unwanted binding sites must be blocked. In this study a variety of reagents and methods were tested for each individual manual IHC protocol. A base set of reliable solutions and techniques for testing antibodies was identified in this study. These included 3%  $H_2O_2$ , low pH, citrate based HIER reagents in a double-boiler, 5% Immunopure<sup>®</sup> Goat Serum in PBS, Leica's IHC Diluent, antibody incubation for one hour at 37 °C, and the Novolink<sup>™</sup> Polymer Detection System. However, with meticulous tailoring of each antibody protocol, our results contained a variety of reagents that were ultimately selected for each staining procedure.

Identification of equine-reactive antibodies can be challenging due to limited knowledge of commercial antibodies' reactivity with veterinary tissues and lack of development of multiple species-specific reagents. Evaluation of cross-species reactivity of commercially available antibodies, particularly CD antigens, has been largely performed in equine tissues prepared for whole-cell analysis like flow-cytometry (*Johne et al., 1997; Mérand et al., 2003; Terio et al., 2003; Kunisch et al., 2004; Ibrahim et al., 2007*) or on fresh or frozen specimens (*Bilzer et al., 1995; Zeng et al., 1996; Lemos et al., 2008; Härtig et al., 2009*). Large screenings of non-equine derived antibodies have often resulted in limited identification of equine reactive reagents (*Ibrahim et al., 2007; Schnabel et al., 2013; Szabo & Gulya, 2013*). Of the 26 antibodies in this study, six antibodies successfully reacted with FFPE equine tissues. All six antibodies were derived from non-equine antigens. Their reactivity to equine tissue is likely due to conserved epitope targets and the intracellular location of the antigen peptides (*Jones et al., 1993; Ahmed et al., 2007*). If an antigen's peptide sequence is highly conserved across multiple species, the production of antibodies from the immunized host may be inhibited due to immune tolerance. However, because these antigens have intracellular origins, the host (i.e. rabbit or mouse) is more likely to elicit an immune response upon antigen exposure. In this way, these

antibodies against human or swine antigen may successfully be produced and react with equine or other species' tissues.

Commercial macrophage antibodies are non-specific to macrophages and often cross react with monocytes, granulocytes, dendritic cells, and fibroblasts (*Johne et al., 1997; Kunisch et al., 2004; Shaw et al., 2005; Sellner et al., 2014*). Multiple macrophage-directed antibodies with different target antigens were tested so that a distinction could be made between tissue macrophage populations and cells that may cross-react with those antibodies. RAM11 (Dako, Glostrup, Denmark) recognizes an uncharacterized, cytoplasmic antigen specific to rabbit macrophages. AM-3K (TransGenic, Strasborg, France) anti-macrophage antibody raised against human alveolar macrophage antigen recognizes cytoplasmic membrane epitopes. Anti-CD68<sup>+</sup> (KP1; Leica, Wetzlar, Germany) recognizes primarily lysosomal membrane proteins of macrophages, secondarily macrophage membranes, and is found on monocytes, neutrophils, basophils, and large lymphocytes. None of these three antibodies were successfully reactive in FFPE brain or lymphoid tissues. Although frequently used as a macrophage-characterizing marker, MAC387<sup>+</sup> (Leica, Wetzlar, Germany) recognizes cytoplasmic, human leukocyte antigen, L1 protein found on neutrophils, monocytes, and certain reactive macrophages. Only MAC387<sup>+</sup> antibody was reactive with FFPE equine tissues, so the positively stained population must be considered as a mix of both macrophages and neutrophils. This cell population had a distinct distribution in serial sections of WNV-infected brain each immunolabeled for MAC387<sup>+</sup>, CD3<sup>+</sup> T lymphocytes, and microglia. In *S. neurona* infected equine tissues, where multi-nucleated giant cells are observed, MAC387<sup>+</sup> antibody did not react with these cells even though there was reaction with single cells within each lesion.

Microglia, which are of a monocytic lineage, were identified with anti-Iba-1-antibody (Wako, Neuss, Germany). This antibody recognizes an intracellular, calcium binding protein basally expressed by both microglia and macrophages, and upregulated in microglia following injury. Cell morphology, specifically microglial processes, may be used to manually separate microglia and macrophage populations. A limited number of equine studies have investigated microglial activation using major histocompatibility complex (MHC) markers (*Mullen, Buck & Smith, 1992; Lemos et al., 2008*). MHC predominantly targets activated microglia populations, whereas anti-Iba-1<sup>+</sup> antibody will stain microglia in both their resting and activated form (*Szabo & Gulya, 2013*). Additionally, using MHC as a marker for microglial cells in an inflammatory disease would also identify all peripheral and CNS cells that have MHC upregulated due to this biological process. Characterization of the total, both resting and activated, microglia population is therefore improved with the use of Iba-1<sup>+</sup> antibody (*Ahmed et al., 2007*).

Regarding detection of neurons, it must be recognized that NF-H proteins are concentrated within axons and are weakly present within cell bodies. Changes in staining intensity within cell bodies and dendrites with NF-H antibodies, such as NAP4 that targets phosphorylated NF-H proteins, are useful for identifying neuropathology (*Shaw et al., 2005; Sellner et al., 2014*). For quantification of neuronal populations, an antibody specific to neuronal nuclei should be used. Neu-N antibodies, which target DNA-binding proteins

of neural nuclei, are typically used for the staining and counting of neuron perikaryons (*Mullen, Buck & Smith, 1992*).

This study successfully identified a collection of protocols for the characterization of cell populations in the histopathology of encephalitic diseases in equine CNS tissues. This suite included IHC staining of CD3<sup>+</sup> T lymphocytes, CD79 $\alpha$ cy<sup>+</sup> B lymphocytes, MAC387<sup>+</sup> macrophages and neutrophils, Iba-1<sup>+</sup> microglia, GFAP<sup>+</sup> astrocytes, and NF-H<sup>+</sup> neurons in FFPE tissues. The suite was expanded with the inclusion of CD4<sup>+</sup> and CD8<sup>+</sup> T lymphocytes in fresh, frozen tissue sections. Useful application of these antibodies was supported with the characterization of neuropathology in diseased horses with clinical and experimental West Nile encephalitis, clinical equine protozoal myeloencephalitis, and, preliminarily, clinical Eastern equine encephalitis. Western blot analysis utilizing these antibodies should next be performed to support the specificity of their reactivity in the horse.

## ACKNOWLEDGEMENTS

We would like to thank Dr. Gerry Shaw and Dr. W. Jake Streit for their consultation and loan of antibodies.

## ADDITIONAL INFORMATION AND DECLARATIONS

### Funding

This work was funded by the Fern Audette Endowment for Equine Studies. The funders had no role in study design, data collection and analysis, decision to publish, or preparation of the manuscript.

### Competing Interests

The authors declare that they have no competing interests.

### Author Contributions

- Gretchen H. Delcambre conceived and designed the experiments, performed the experiments, analyzed the data, wrote the paper, prepared figures and/or tables, reviewed drafts of the paper.
- Junjie Liu performed the experiments, analyzed the data, reviewed drafts of the paper.
- Jenna M. Herrington performed the experiments.
- Kelsey Vallario performed the experiments.
- Maureen T. Long conceived and designed the experiments, performed the experiments, analyzed the data, contributed reagents/materials/analysis tools, wrote the paper, reviewed drafts of the paper.

### Animal Ethics

The following information was supplied relating to ethical approvals (i.e., approving body and any reference numbers):

Archived equine tissues were utilized under the University of Florida (Gainesville, FL) Institutional Animal Care and Use Committee protocols #F077, #F093, #D163, and #4109.

## Data Deposition

The following information was supplied regarding data availability:

The research in this article did not generate any raw data.

## Supplemental Information

Supplemental information for this article can be found online at <http://dx.doi.org/10.7717/peerj.1601#supplemental-information>.

## REFERENCES

- Ahmed Z, Shaw G, Sharma VP, Yang C, McGowan E, Dickson DW. 2007.** Actin-binding proteins coronin-1a and IBA-1 are effective microglial markers for immunohistochemistry. *Journal of Histochemistry and Cytochemistry* **55**(7):687–700 DOI [10.1369/jhc.6A7156.2007](https://doi.org/10.1369/jhc.6A7156.2007).
- Beckstead JH. 1994.** A simple technique for preservation of fixation-sensitive antigens in paraffin-embedded tissues. *Journal of Histochemistry and Cytochemistry* **42**(8):1127–1134 DOI [10.1177/42.8.8027531](https://doi.org/10.1177/42.8.8027531).
- Bilzer T, Planz O, Lipkin WI, Stitz L. 1995.** Presence of CD4<sup>+</sup> and CD8<sup>+</sup> T cells and expression of MHC class I and MHC class II antigen in horses with Borna disease virus-induced encephalitis. *Brain Pathology* **5**(3):223–230 DOI [10.1111/j.1750-3639.1995.tb00598.x](https://doi.org/10.1111/j.1750-3639.1995.tb00598.x).
- Bourgeois MA, Denslow ND, Seino KS, Barber DS, Long MT. 2011.** Gene expression analysis in the thalamus and cerebrum of horses experimentally infected with West Nile virus. *PLoS ONE* **6**(10):e24371 DOI [10.1371/journal.pone.0024371](https://doi.org/10.1371/journal.pone.0024371).
- Boy MG, Galligan DT, Divers TJ. 1990.** Protozoal encephalomyelitis in horses: 82 cases 1972–1986. *Journal of the American Veterinary Medical Association* **196**(4):632–634.
- Buchwalow I, SamoiloVA V, Boecker W, Tiemann M. 2011.** Non-specific binding of antibodies in immunohistochemistry: fallacies and facts. *Scientific Reports* **1**:28 DOI [10.1038/srep00028](https://doi.org/10.1038/srep00028).
- Cantile C, Del Piero F, Di Guardo G, Arispici M. 2001.** Pathologic and immunohistochemical findings in naturally occurring West Nile virus infection in horses. *Veterinary Pathology* **38**(4):414–431 DOI [10.1354/vp.38-4-414](https://doi.org/10.1354/vp.38-4-414).
- Daneshtalab N, Doré JJE, Smeda JS. 2010.** Troubleshooting tissue specificity and antibody selection: procedures in immunohistochemical studies. *Journal of Pharmacological and Toxicological Methods* **61**(2):127–135 DOI [10.1016/j.vascn.2009.12.002](https://doi.org/10.1016/j.vascn.2009.12.002).
- Dubey JP, Lindsay DS, Saville WJ, Reed SM, Granstrom DE, Speer CA. 2001.** A review of Sarcocystis neurona and equine protozoal myeloencephalitis (EPM). *Veterinary Parasitology* **95**(2–4):89–131 DOI [10.1016/S0304-4017\(00\)00384-8](https://doi.org/10.1016/S0304-4017(00)00384-8).
- Ferrier CM, van Geloof WL, de Witte HH, Kramer MD, Ruiter DJ, van Muijzen GN. 1998.** Epitopes of components of the plasminogen activation system are re-exposed in formalin-fixed paraffin sections by different retrieval techniques. *Journal of Histochemistry and Cytochemistry* **46**(4):469–476 DOI [10.1177/002215549804600406](https://doi.org/10.1177/002215549804600406).
- Grosche A, Freeman DE, Morton AJ, Polyak MMR, Matyjaszek SA. 2012.** Effects of ischemia and reperfusion on production of nitrotyrosine, activation of eosinophils, and apoptosis in the large colonic mucosa of horses. *American Journal of Veterinary Research* **73**(1):53–61 DOI [10.2460/ajvr.73.1.53](https://doi.org/10.2460/ajvr.73.1.53).
- wGutierrez M, Forster FI, McConnell SA, Cassidy JP, Pollock JM, Bryson DG. 1999.** The detection of CD2<sup>+</sup>, CD4<sup>+</sup>, CD8<sup>+</sup>, and WC1<sup>+</sup> T lymphocytes, B cells and macrophages in fixed and paraffin embedded bovine tissue using a range of antigen recovery and signal amplification techniques. *Veterinary Immunology and Immunopathology* **71**(3–4):321–334 DOI [10.1016/S0165-2427\(99\)00106-3](https://doi.org/10.1016/S0165-2427(99)00106-3).

- Härtig W, Reichenbach A, Voigt C, Boltze J, Bulavina L, Schuhmann MU, Seeger J, Schusser GF, Freytag C, Grosche J. 2009. Triple fluorescence labelling of neuronal, glial and vascular markers revealing pathological alterations in various animal models. *Journal of Chemical Neuroanatomy* 37(2):128–138 DOI 10.1016/j.jchemneu.2008.10.003.
- Ibrahim S, Saunders K, Kydd JH, Lunn DP, Steinbach F. 2007. Screening of anti-human leukocyte monoclonal antibodies for reactivity with equine leukocytes. *Veterinary Immunology and Immunopathology* 119(1–2):63–80 DOI 10.1016/j.vetimm.2007.06.034.
- Johne B, Fagerhol MK, Lyberg T, Prydz H, Brandtzaeg P, Naess-Andresen CF, Dale I. 1997. Functional and clinical aspects of the myelomonocyte protein calprotectin. *Molecular Pathology* 50(3):113–123 DOI 10.1136/mp.50.3.113.
- Jones M, Cordell JL, Beyers AD, Tse AG, Mason DY. 1993. Detection of T and B cells in many animal species using cross-reactive anti-peptide antibodies. *Journal of Immunology* 150(12):5429–5435.
- Krenacs L, Krenacs T, Stelkovics E, Raffeld M. 2010. Heat-induced antigen retrieval for immunohistochemical reactions in routinely processed paraffin sections. *Methods in Molecular Biology* 588:103–119 DOI 10.1007/978-1-59745-324-0\_14.
- Kunisch E, Fuhrmann R, Roth A, Winter R, Lungershausen W, Kinne RW. 2004. Macrophage specificity of three anti-CD68 monoclonal antibodies (KP1, EBM11, and PGM1) widely used for immunohistochemistry and flow cytometry. *Annals of the Rheumatic Diseases* 63(7):774–784 DOI 10.1136/ard.2003.013029.
- Kumar GL, Rudbeck L. 2009. Chapter 8 Demasking of antigens. In: Kumar GL, Rudbeck L, eds. *Immunohistochemical Staining Methods*. Fifth edition. Carpinteria: Dako North America, 51–58.
- Lemos KR, Marques LC, Aquino LPCT, Alessi AC, Zacarias RZ. 2008. Astrocytic and microglial response and histopathological changes in the brain of horses with experimental chronic *Trypanosoma evansi* infection. *Revista do Instituto de Medicina Tropical de São Paulo* 50(4):243–249 DOI 10.1590/s0036-46652008000400011.
- Matyjaszek SA, Morton AJ, Freeman DE, Grosche A, Polyak MMR, Kuck H. 2009. Effects of flunixin meglumine on recovery of colonic mucosa from ischemia in horses. *American Journal of Veterinary Research* 70(2):236–246 DOI 10.2460/ajvr.70.2.236.
- Mérant C, Bonnefont C, Desbos A, Greenland T, Cadore J-L, Monier J-C. 2003. Cross-species reactivity of seven monoclonal antibodies with equine lymphocytes by flow cytometry. *Veterinary Research* 34(6):791–801 DOI 10.1051/vetres:2003033.
- Mullen RJ, Buck CR, Smith AM. 1992. NeuN, a neuronal specific nuclear protein in vertebrates. *Development* 116(1):201–211.
- Porter MB, Long MT, Getman LM, Giguère S, MacKay RJ, Lester GD, Alleman AR, Wamsley HL, Franklin RP, Jacks S, Buergelt CD, Detrisac CJ. 2003. West Nile virus encephalomyelitis in horses: 46 cases (2001). *Journal of the American Veterinary Medical Association* 222(9):1241–1247 DOI 10.2460/javma.2003.222.1241.
- Ramos-Vara JA. 2005. Technical aspects of immunohistochemistry. *Veterinary Pathology* 42(4):405–426 DOI 10.1354/vp.42-4-405.
- Ritchey JW, Black DH, Rogers KM, Eberle R. 2006. In vivo experimentation with simian herpesviruses: assessment of biosafety and molecular contamination. *Journal of the American Association for Laboratory Animal Science* 45(2):7–12.
- Rushton JO, Lecollinet S, Hubálek Z, Svobodová P, Lussy H, Nowotny N. 2013. Tick-borne encephalitis virus in horses, Austria, 2011. *Emerging Infectious Diseases* 19(4):635–637 DOI 10.3201/eid1904.121450.

- Schnabel CL, Wagner S, Wagner B, Durán MC, Babasyan S, Nolte I, Pfarrer C, Feige K, Murua Escobar H, Cavalleri J-MV. 2013. Evaluation of the reactivity of commercially available monoclonal antibodies with equine cytokines. *Veterinary Immunology and Immunopathology* 156(1–2):1–19 DOI 10.1016/j.vetimm.2013.09.012.
- Seino KK, Long MT, Gibbs EPJ, Bowen RA, Beachboard SE, Humphrey PP, Dixon MA, Bourgeois MA. 2007. Comparative efficacies of three commercially available vaccines against West Nile Virus (WNV) in a short-duration challenge trial involving an equine WNV encephalitis model. *Clinical and Vaccine Immunology* 14(11):1465–1471 DOI 10.1128/CVI.00249-07.
- Sellner J, Davies NW, Howard RS, Petzold A. 2014. Neurofilament heavy chain as a marker of neuroaxonal pathology and prognosis in acute encephalitis. *European Journal of Neurology* 21(6):845–850 DOI 10.1111/ene.12390.
- Shaw G, Yang C, Ellis R, Anderson K, Parker Mickle J, Scheff S, Pike B, Anderson DK, Howland DR. 2005. Hyperphosphorylated neurofilament NF-H is a serum biomarker of axonal injury. *Biochemical and Biophysical Research Communications* 336(4):1268–1277 DOI 10.1016/j.bbrc.2005.08.252.
- Shi SR, Key ME, Kalra KL. 1991. Antigen retrieval in formalin-fixed, paraffin-embedded tissues: an enhancement method for immunohistochemical staining based on microwave oven heating of tissue sections. *Journal of Histochemistry and Cytochemistry* 39(6):741–748 DOI 10.1177/39.6.1709656.
- Szabo M, Gulya K. 2013. Development of the microglial phenotype in culture. *Neuroscience* 241:280–295 DOI 10.1016/j.neuroscience.2013.03.033.
- Terio KA, Stalis IH, Allen JL, Stott JL, Worley MB. 2003. Coccidioidomycosis in Przewalski's horses (*Equus przewalskii*). *Journal of Zoo and Wildlife Medicine* 34(4):339–345 DOI 10.1638/02-005.
- van Marle GG, Antony JJ, Ostermann HH, Dunham CC, Hunt TT, Halliday WW, Maingat FF, Urbanowski MDM, Hobman TT, Peeling JJ, Power CC. 2007. West Nile virus-induced neuroinflammation: glial infection and capsid protein-mediated neurovirulence. *Journal of Virology* 81(20):10933–10949 DOI 10.1128/JVI.02422-06.
- Wendelboe HG, Bisgaard K. 2013. Chapter 17-Background. In: Kumar GL, Rudbeck L, eds. *Immunohistochemical Staining Methods*. Fifth edition. Carpinteria: Dako. Available at [http://www.dako.com/08002\\_ihc\\_staining\\_methods.pdf](http://www.dako.com/08002_ihc_staining_methods.pdf).
- Yu G-Y, Wiley MR, Kugelman JR, Ladner JT, Beitzel BF, Eccleston LT, Morazzani EM, Glass PJ, Palacios GF. 2015. Complete coding sequences of eastern equine encephalitis virus and venezuelan equine encephalitis virus strains isolated from human cases. *Genome Announcements* 3(2):e00243–15 DOI 10.1128/genomeA.00243-15.
- Zeng L, Takeya M, Ling X, Nagasaki A, Takahashi K. 1996. Interspecies reactivities of anti-human macrophage monoclonal antibodies to various animal species. *Journal of Histochemistry and Cytochemistry* 44(8):845–853 DOI 10.1177/44.8.8756757.